

Phone +91 22 42184218 Fax +91 22 42184350

Email tcmktg@sarex.com URL www.sarex.com

## PRODUCT LITERATURE

### SARAGUARD-AM (CONC)

#### Product Information

Saraguard-AM (Conc) is a concentrated antimicrobial agent for finishing of textile substrate. Saraguard-AM (Conc) can be applied by padding, exhaust soaking and spraying methods. Saraguard-AM (Conc) is non-leaching antimicrobial agent. It is effective against a broad spectrum of microbes. It shows antibacterial and antifungal activity by contact mechanism.

It can be applied by spray technique on the surface of shoes and packing box or paper to prevent mildew formation.

#### Key Features & Benefits

Key features  Non-leaching	Benefits  Durable to multiple home launderings
Concentrated	Economical
Broad spectrum	Works on all types of microbes such as Staphylococcus aureus, Escherichia coli, Klebsiella pneumoniae, Methicillin- Resistant Staphylococcus aureus (MRSA).
Universal	Suitable for cotton, polyester, polyamide and acrylic

#### General Characteristics

Physical appearance	:	Yellow to Orange, transparent to hazy liquid
Ionic nature	:	Cationic
pH of 1% solution	:	5 +/- 1
Miscibility	:	Miscible with water
Compatibility	:	Compatible with cationic and nonionic
		products
Stability	:	Stable to dilute acids and unstable to dilute alkalies



Phone +91 22 42184218 Fax +91 22 42184350

Email tcmktg@sarex.com URL www.sarex.com

# PRODUCT LITERATURE

### SARAGUARD-AM (CONC)

#### Application

Padding proces	S			
Saraguard-AM ( Saralink-ULF Pick-up Bath pH Drying	:	<pre>10-40 g/l 30 g/l + 6 g/l MgC12 (Optional,to improve extended durability) 65-70 % 5.0-6.0 140-160 deg.C, 1 min. (all substrates)</pre>		
Exhaust proces				
Saraguard-AM ( Bath pH Bath Temp. Time	Conc)* : :	1-4 % 5.0-6.0 30-40 deg.C 15-20 mins.		
* Dosages reco strength.	mmended a	are based on the concentrated products		
Dosages of d solid conten	_	coduct to be calculated on the basis of		
Precautions:				
The product is supplied in the anhydrous form. The product to be diluted with water just before a treatment only when ready for processing of textiles. Once added to the pad bath / exhaust bath the stability of such a bath is limited; due to hydrolysis of the active ingredients in our antimicrobial leading to self cross linking. Hence, padding bath containing Saraguard-AM (Conc) should be used up in the shortest time possible, as extended storage of such bath would lead to formulation of insoluble particles which may settle as precipitation and reduces performance.				
Kindly note that the product is readily dilutable and stored in anhydrous polar solvents like methanol, ethanol, iso-propyl alcohol or low molecular weight glycol like MEG or DEG.				

Test Results:



Phone +91 22 42184218 Fax +91 22 42184350

Email tcmktg@sarex.com URL www.sarex.com

## PRODUCT LITERATURE

## SARAGUARD-AM (CONC)

- Cotton, polyester and their blended fabrics finished with 30 g/l Saraguard-AM (Conc) shows 80-95% reduction in bacterial count by AATCC-100, after 20 Home Launderings.
- Polyamide and Acrylic fabrics finished with 30 g/l Saraguard-AM (Conc) shows 95-100% reduction in bacterial count by AATCC 100 method.

#### Instructions For Dilution

Procedure to prepare 40-48% diluted product of 70% Saraguard-AM (Conc)

Dilution	40%	48%
Saraguard-AM (Conc)	57	68
Dry Butyl Di-glycol OR Dry	43	32
Methanol or Dry Ethanol		
Total	100	100

- Take required quantity of Saraguard-AM (Conc)
- To this, gradually add any one of the anhydrous solvent mentioned above under constant stirring to obtain homogeneous product

Note:

Even a drop of water will polymerise the material. Care should be taken to avoid water contact at any stage of dilution or during packing.

#### Precautions

Storage	: Store in cool, ventilated shed away from
	heat and direct sunlight.Storage temperature should not exceed 35 deg C.Close lids firmly to avoid contact with air and moisture.
Shelf Life	: 9 months from the date of manufacturing, if stored under controlled conditions.

[The above information is given in good faith and is without warranty] [Prod.Code :PY731 | Last Upd. On : 04/11/15 ] [PP]



Phone +91 22 42184218 Fax +91 22 42184350

Email tcmktg@sarex.com URL www.sarex.com

# PRODUCT LITERATURE

SARAGUARD-AM (CONC)

This is a computer generated report. Hence not signed.