

501 Waterford, C Wing, C.D.Barfiwala Marg(Juhu Lane), Andheri(W),Mumbai 400 058. INDIA

Phone +91 22 42184218 Fax +91 22 42184350

Email tcmktg@sarex.com URLwww.sarex.com

PRODUCT LITERATURE

HELOPREP-GC

Product Information

Heloprep-GC is a novel biodegradable all-in-one product which is suitable for bleaching of cellulosic substrates and its blends. Heloprep-GC possess sequestering and peroxide stabilization properties, hence addition of sequestering agent and peroxide stabilizer is not required in treatment bath. Bleaching with Heloprep-GC exhibits very good softness, absorbency, rewetting and whiteness of bleached fabric than conventional bleached fabric. Yarn bleached with Heloprep-GC decreases the coefficient of friction, improves strength and reduces hariness to a certain extent. Suitable for Autodosing / dispensing.

Key Features & Benefits

Key Features	Benefits
Low alkali during bleaching	No separate rinsing after bleaching. Saving in water
In built sequestering peroxide stabilization	No need to add sequestering agent or peroxide stabilizer
Eco-friendly	Low BOD/COD and TDS in effluent bath

General Characteristics

Physical appearance	:	Light Blue liquid
Ionic nature	:	Nonionic to weakly anionic
pH of 1% solution	:	3.0 to 3.5
Miscibility	:	Misible with water
Compatibility	:	Compatible with anionic and non-ionic products
Stability	:	Stable to dilute acids and dilute alkalies

Application

Recipe for RFD (Ready for	dyeing)
Heloprep-GC	: 0.8-1.0%
Saracrease-HG(Conc)	: 0.1%
Caustic flakes	: 1.5-1.8%
Hydrogen Peroxide(50% v/v)	: 2.0-2.5%
Temperature	: 95-98 deg.C for 30-45min.
	Heloprep-GC Saracrease-HG(Conc) Caustic flakes Hydrogen Peroxide(50% v/v)



501 Waterford, C Wing, C.D.Barfiwala Marg(Juhu Lane), Andheri(W),Mumbai 400 058. INDIA

Phone +91 22 42184218 Fax +91 22 42184350

Email tcmktg@sarex.com URLwww.sarex.com

PRODUCT LITERATURE

HELOPREP-GC

2. Recipe for full white Heloprep-GC : 1.2-1.5% Saracrease-HG (Conc) : 0.1% Caustic flakes : 2.25 - 2.75% Hydrogen Peroxide : 6.0-8.0% Temperature : 95-98 deg.C for 60-75min.

Optical brightener to be dosed in hot wash OR in the same bleached bath after 50-60 min.

Precautions

Storage	: Stor	e in c	cool, v	ventilat	ed she	ed away	from
				2	-	ge temper	
	shou	ld not	exceed	35 deg	C.Clos	se lids f	irmly
	to	avoid	contact	t with	air	and mois	ture.

Shelf Life : 9 months from the date of manufacturing, if stored under controlled conditions.

[The above information is given in good faith and is without warranty] [Prod.Code :L009226 | Last Upd. On : 11/03/24] [PVR]

This is a computer generated report. Hence not signed.